

OBSTETRICS

innovate. precision.

***Birth
Assistance***



!theronx

Quality Standard

All Medical instruments are manufactured under strict quality control procedures to ensure consistent performance, safety, and reliability. Each instrument is produced using high-grade medical stainless steel and undergoes precise machining, heat treatment, and finishing to achieve optimal strength, sharpness, and corrosion resistance.

Clean Precision

Every production stage is carefully monitored, from raw material selection to final inspection. Instruments are tested for accuracy, durability, smooth operation, and surface finish to meet international medical standards. Special attention is given to edge retention, joint alignment, and ergonomic balance to ensure dependable performance during surgical procedures.

ENGINEERED FOR
CLEAN PRECISIONS



15A	15A	15A	25A	15A
Theronx.	www.theronx.com	+92 327 911 7077	!theronx	info@theronx.com



Global Reach & Customer Satisfaction

Global Reach

Theronx has a global presence. We serve customers from various countries, ensuring that medical professionals and beauty experts worldwide have access to premium surgical and Medical instruments.

Trustworthiness

We value the trust our customers place in us and continuously strive to exceed their expectations. Your feedback is valuable to us, and we use it to improve our products and services to serve you better.

Customer Satisfaction

At Theronx “Customer Satisfaction” is our top priority. We are committed to providing an exceptional shopping experience and unparalleled customer support. Our knowledgeable and friendly team is always available to assist you, answer your queries, and guide you in selecting suitable instruments for your specific requirements.

ABOUT us

Crafting Precision.

Theronx is a Pakistan-based manufacturer and supplier of premium surgical and Medical instruments, serving global medical and surgical markets from Sialkot, the world-renowned hub of surgical craftsmanship.

Building Trust.

Since 2010, we have been dedicated to producing instruments that combine precision engineering, durable materials, and consistent quality, meeting the expectations of medical professionals and healthcare experts worldwide.

We take pride in being a Pakistani company, contributing to the legacy of Sialkot's surgical industry while upholding international manufacturing standards.

A Legacy of Quality & Innovation

With over a decade of industry experience, Theronx has established itself as a dependable partner for healthcare and medical brands. Our manufacturing process is driven by skilled craftsmanship, modern machinery, and strict quality control ensuring every instrument reflects our commitment to excellence.

By understanding the real-world needs of professionals, we deliver tools that perform reliably in demanding environments.





Contact Info

Email

info@theronx.com
sales@theronx.com

Phone

+92-327-911-7077

Office Address:

Theronx INC,
Muhammad,
pura, Sialkot

Pelvimeters

Pelvimeter
Pelvimetros
Pelvimeters



1/2

Martin
21-100-35
35 cm
cm / inch scale



1/2

Collin
21-102-33
33 cm / inch scale



1/2

Collin
21-104-35
35 cm / inch scale

The additional straight arm can be replaced. The pelvimeter with straight arm is used for localisation of the portio vaginalis uteri and of vaginal metastasis.





Breisky
21-106-33
33 cm
Graduated in centimeters



De Lee
21-108-20
20 cm / inch scale

For internal measurements
between symphysis and
promontory

Obstetrical Forceps

Geburtszange

Pinça Obstétrica

Pinzas obstétricas



21-130-36

36.5 cm

21-130-40

40.5 cm

Nägele

21-130-36 - 21-130-40

36.5 cm



Wrigley

21-132-28

28 cm



Simpson

21-134-30

30.5 cm




Simpson

21-136-36

36.5 cm



Piper
21-138-44
44.5 cm



Kielland
21-140-41
41 cm



Kielland-Luikart
21-142-40
40 cm

Obstetrical Forceps

Geburtszange

Pinça Obstétrica

Pinzas obstétricas



Tarnier
21-144-40
40 cm




Braun
21-146-42
42 cm
Cranioclast

Obstetrical Hooks, Perforator
Geburtshelfer, Perforator
Ganchos Obstétricos, Perforador
Ganchos Obstétricos, Perforador



Smellie
21-170-30
30.5 cm
Obstetrical hook



Braun
21-172-31
31 cm
Obstetrical hook




Smellie
21-178-27
27 cm
Perforator

Placenta, Ovum Forceps

Plazenta, Ovum Pinzette

Placenta, Fórceps Ovum

Placenta, fórceps de óvulo



Winter
21-190-01 - 21-190-03
29 cm
Straight




Fig: 1
21-190-01



Fig: 2
21-190-02



Fig: 3
21-190-03



Winter
21-192-01 - 21-192-03
29 cm
Curved




Fig: 1
21-192-01



Fig: 2
21-192-02



Fig: 3
21-192-03




Haemorrhoidal Forceps, Polypus, Ovum Forceps

Hämorrhoiden, Polypus, Ovum Pinzette

Pinça Hemorroidal, Polypus, Pinça Ovum

Pinzas hemorroidales, pólipos, pinzas de óvulos




Luer
21-206-22 - 21-206-24
22 cm
Haemorrhoidal forcep

Doyen
21-208-18
18.5 cm
Polypus and ovum forcep


Bumm
21-210-28
28 cm
Polypus and ovum forcep

Placenta, Ovum Forceps


Plazenta, Ovum Pinzette
Placenta, Fórceps Ovum
Placenta, fórceps de óvulo




Green-Armytage
21-220-22
21 cm
Uterine clamp



Sanger
21-222-27 - 21-224-27
27.5 cm
Placenta and ovum forcep



21-222-27
27.5 cm
Straight



21-224-27
27 cm
Curved

Scalp Flap Forceps, Obstetrical Lever

Kopfhaut-Flap-Zange, Geburtshelfer

Pinça Flap do couro cabeludo, alavanca obstétrica


Pinzas de aleta de cuero cabelludo, palanca obstétrica




Willett
21-236-18
18.5 cm
Scalp flap forceps

Gauss
21-238-26
26 cm
Scalp flap forceps


Sellheim
21-246-32
32 cm
Obstetrical lever




Back side serrated



21-254-01
14,5 mm



21-254-02
18 mm



21-254-03
20 mm

Cuzzi
21-254-01 - 21-254-03
30 cm
Blunt

Placenta Scoops, Umbilical Cord Clamps


Placenta Schaufeln, Nabelschnurklemmen

Colheres de placenta, braçadeiras de cordão umbilical

Placenta Scoops, abrazaderas del cordón umbilical



Pestalloza
21-256-30
30.5 cm
Placenta scoop



Kane
21-264-08
8.5 cm
Umbilical cord clamp



Collin
21-266-08
8.5 cm
Umbilical cord clamp